

## E.max ZirCAD Colouring Liquid

Staining solution for IPS e.max ZirCAD MT by way of brush infiltration before the sintering process. These colouring liquids are available in 16 A-D tooth shades for easy and individual customisation.

## E.max ZirCAD Colouring Liquid Indicator

Take the guess work out of zirconia stain infiltration with these liquids. ZirCAD colour Liquid concentrates can be mixed with ZirCAD MT and LT colouring liquids to render the infiltration of ZirCAD restorations visible.

Shade	15ml
Yellow	JDI852
Red	JDI853
Blue	JDI854

## E.max ZirCAD Colouring Liquid Diluter

The IPS e.max ZirCAD Colouring Liquid Diluter is a liquid for diluting the IPS e.max ZirCAD Colouring Liquids. This allows colouring liquids with a lower shading power to be produced, which results in lighter tones.

60ml
JDI855



Size	A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4	Blue	Violet	Grey	Orange	Brown	Purple
	MT (Medium Translucency)																					
60ml	JDI831	JDI832	JDI833	JDI834	JDI835	JDI836	JDI837	JDI838	JDI839	JDI840	JDI841	JDI842	JDI843	JDI844	JDI845	JDI846	-	-	-	-	-	-
15ml	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	JDI847	JDI848	JDI849	JDI850	JDI851	-
	LT (Low Translucency)																					
60ml	JDI809	JDI810	JDI811	JDI812	JDI813	JDI814	JDI815	JDI816	JDI817	JDI818	JDI819	JDI820	JDI821	JDI822	JDI823	JDI824	-	-	-	-	-	-
15ml	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	-	JDI825	JDI826	JDI827	JDI828	JDI829	JDI830