

E.max CAD Shades

Dentin stains in paste form that are optimally coordinated with IPS e.max CAD. Can be applied both on “blue” or crystallized restorations.

Shade

3g

1	JDI221
2	JDI222
3	JDI223
4	JDI224
11	JDI225
12	JDI226



E.max CAD Stains

Intensive stains in paste form. They may be applied both on “blue” or crystallized restorations.

Shade

1g

Creme	JDI228
Khaki	JDI232
Mahagoney	JDI233
Olive	JDI231
White	JDI227
Sunset	JDI229
Copper	JDI230



E.max CAD Glaze Liquid 15ml

Special liquid to dilute IPS e.max CAD Crystall./Shades and Stains and the Glaze Paste.

JDI237



E.max CAD Glaze Paste 3g

Special glaze paste to be used with fully anatomical IPS e.max CAD restorations. It may be applied either prior to crystallization on the "blue" crown or in the subsequent firing cycle.

JDI234



E.max CAD Glaze Spray 270ml

Special glaze spray to be applied on restorations once IPS e.max CAD Crystall Shades and Stains have been applied. The restorations are fired in a combined crystallisation and glaze firing cycle.

JDI235



E.max CAD Add-On - Incisal 5g

For adjustments, such as proximal or occlusal contact points. These adjustments can be made directly with the crystallization firing or in a subsequent corrective firing cycle.

JDI236



E.max CAD Crystallization Tray

Silicon nitride firing tray for the crystallization of IPS e.max CAD restorations. The Crystallization pins supplied are required for the crystallization of IPS e.max CAD LT restorations.

JDI239



E.max CAD Crystallization Pins

Refill of Small, Medium and Large pins for use with Crystallization Tray.

JDI240

